

FOLD & GO

COMPACT FOLDING POWERCHAIR
P113



 **merits**

AUSTRALIA

www.meritsaustralia.com












FEATURES

- Lightweight, foldable power wheel chair which folds easily in seconds for transportation or storage
- Folds in one easy step
- Compact size means great manoeuvrability in tight spaces
- Removable battery box with battery charge indicator & extra charge port for charging off board
- Adjustable height footplate
- Adjustable height backrest
- Intelligent, regenerative, electromagnetic braking
- Vibrant smart colours available, red, blue or grey

SPECIFICATION

Max. Range	11km / 7mi. 15km / 9.6mi. (optional)
Max. Speed	7.5kph / 4.7kph
Battery	1pc Lithium-ion battery 24V/11.4AH / 17.1AH (optional)
Motor	DC 24V, 70W (Rating)/210W (Max)
Gradient	3°
Controller	PG VR2 50A
Charger	2A
Caster	7" PU tire
Drive Wheel	8"(200x50) PU tire
Brake	electromagnetic brakes

DIMENSIONS (cm/in)

 Seat Width	 Seat Depth	 Seat Height from Groud	
46/18	42/16.5	52/20.5	
 Armrest Height	 Backrest Height	 Footrest Extension	 Overall Length
22/9	33/13	42/16.5	48/19
 Overall Length	 Overall Width	 Overall Height	 Folded Dimension
97cm/38"	59/23	95-100/ 37.5-39	59x37x81/ 23x15x32 LxWxH

WEIGHT (kg/lb)

ttl w/batt.	ttl w/o batt.	capacity
23/50.5 23.9/52.5 (optional)	20.9/46	114/250

COMPACT FOLDING POWER WHEELCHAIR



 **merits**
AUSTRALIA

www.meritsaustralia.com

Merits Australia is committed to constantly improve and develop products, we reserve the right to alter specifications without prior notice. Pictures for illustration purposes only.

FDA

CE

MDR 2017/755



ISO 13485

ISO 14001

rev2018-1